


# SGTi-flex COVID-19 IgG



**FDA Emergency Use Authorized.**  
Healthcare professional use only

## COVID-19 IgG test?

- ▶ IgG antibody was generated for COVID-19 within 1 week after symptom onset and remains in the blood after recovery. Moreover, IgG has much better neutralizing activity and the so more elongated expression.
- ▶ IgG tests is intended for use as an aid in identifying individuals with an adaptive immune response to SARSCoV-2, indicating recent or prior infection.



Hou H, Wang T, Zhang B et al., Detection of IgM and IgG antibodies in patients with coronavirus disease 2019. Clin Transl Immunology 2020; 9(5): e01136.

## IgG and Immunity

- ▶ IgG tests can help determine the proportion of a population previously infected with SARS-CoV-2 and provide information about populations that may be immune and potentially protected.
- ▶ Serologic test results may assist with identifying persons potentially infected with SARS-CoV-2 and determining who may qualify to donate blood that can be used to manufacture convalescent plasmaexternal icon as a possible treatment for those who are seriously ill from COVID-19.

**FDA Emergency Use Authorized.**



[25 tests/kit]

## performance

### Clinical Performance analysis

TOTAL		Real Time RT-PCR		
		Pos	Neg	Total
SGTi-flex COVID-19 IgG	Pos	171	2	173
	Neg	14	232	246
Total		185	234	419

- Sensitivity : 92.43% (171/185, 95% CI: 87.70%~95.44%)
- Specificity : 99.15% (232/234, 95% CI: 96.94%~99.77%)

### Positive Agreement by Days Post-Symptom Onset

Days from symptom onset	Number of samples tested	IgG Positive Results	IgG PPA (95% CI)
0~7 days	17	7	41.2% (21.61%~63.99%)
8~14 days	24	22	91.7% (74.15%~97.68%)
15~21 days	45	43	95.6% (85.17%~98.77%)
≥ 22 days	99	99	100.00% (96.26%~100.00%)

# SGTi-flex COVID-19 IgG




**FDA Emergency Use Authorized.**  
**Healthcare professional use only**

## Test Procedure

### 1 Collecting of Sample


For the test, 10  $\mu$ l of whole blood, plasma or serum is used.

- Collect fingertip blood using a pipette or blood transfer pipette.
- Or use a blood sample obtained by venipuncture.




### 2 Adding of Sample

Add the collected sample (whole blood/serum/plasma) to the sample well of the test cassette.



### 3 Dropping of Sample buffer

Add 3 drops (90 $\mu$ l) of sample into the sample well of the test cassette.



### 4 Reading Test result

Read the results after 10-15 min.  
The results are invalid after 30 min.




Read after 10 min.




Do not read after 30 min.

## Interpretation


### Positive



### Negative



### Invalid/Retest



1. The test is for in vitro diagnostic use only.
2. As in case of all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single but should rather be made after test all the clinical findings have been evaluated.
3. Negative results do not rule out SARS-CoV-2 infection, particularly in those who have been in contact with the virus. Follow-up testing with a molecular diagnostic should be considered to rule out infection in these individuals.
4. Results from antibody testing should not be used as the sole basis to diagnose or exclude SARS-CoV-2 infection or to inform infection status.
5. Positive results may be due to past or present infection with non-SARS-CoV-2 coronavirus strains, such as coronavirus HKU1, NL63, OC43, or 229E.